

PREMIUM WRISTBANDS

FOR ADULTS AND INFANTS

THE BEST QUALITY & MOST RELIABLE BARCODE WRISTBANDS FOR HEALTHCARE AND ENTERTAINMENT APPLICATIONS

- Economical, durable and water resistant
- Combines ultra comfort with high performance and safety
- Withstands even the harshest hospital environments
- Anti-wrinkle material makes it easy to read barcodes and other identifiers
- Durable non-transferable adhesive closure
- High-quality dark & sharp print that doesn't fade
- Easy to print and apply to patients in seconds
- Water-resistant and tamper-evident



WRISTBANDS

ADULTS

Perfect to withstand the harshest conditions in the hospital environment. Wristbands are made from satin-like material which is soft and gentle, even to the most sensitive skin. Provides a high-quality bar code for scanning and accurate patient ID. Supports text, linear and 2D barcodes, photos and graphics. Features top quality tamper-evident adhesive and antimicrobial additive to protect wristband surface against non-pathogenic bacteria.
Size of Band: 29x292 mm, Information area: 29x114 mm.

APPLICATIONS

HEALTHCARE

- Patient ID tracking
- Hospital admissions
- Hospital nursery

ENTERTAINMENT & HOSPITALITY

- Amusement parks
- Water parks
- Event admissions
- Night club age verification

OTHER MARKETS

- Crowd control
- Animal ID
- Restricted area access

TECHNICAL INFORMATION

	Adult
Product number	WDA-105-OL
Band width	29 mm
Band length	292 mm
Information area	29mm x 114mm
Printing method	Direct Thermal
Material	Satin
Colour	White
Type of closure	Adhesive
Core size	1" / 25 mm
Wound	OUT
Roll O.D.	4,3" / 109mm
Black mark	Yes, at perforation
Bands per roll	150
Packaging	2 rolls/ctn = 300 bands/ctn

PREMIUM WRISTBANDS



Additional features

Phthalate-free, Latex-free, Antimicrobial, Resistant to alcohol, water, soap, hand-sanitizer, and all the most common fluids in the hospital environment.
Certificates: ISO 10993-10 Skin Irritation; Natural Rubber Latex (NRL) free, Bisphenol A (BPA) free, Phthalates (DEHP, DINP, DIDP, DnHP, BBP, DBP, DIBP, DCHP) free.

Wristband safety

Keep wristband clean, dry and don't wear it too tight. When applying a wristband place "one finger under" the band to ensure it is not too tight. When wristbands are worn too tight perspiration, water, soap, lotions etc. caught under the band may cause irritation. If irritation occurs, remove wristband and treat symptoms.

Hand Hygiene: Use proper hand hygiene to reduce the spread of bacteria or other sources of contamination to wristband.

Storage temperature: 22-35°C

TECHNICAL SPECIFICATION



Product specifications are subject to change without notice.