

PREMIUM WRISTBANDS

FOR INFANTS/NEWBORN

THE BEST QUALITY & MOST RELIABLE BARCODE WRISTBANDS FOR HEALTHCARE AND ENTERTAINMENT APPLICATIONS

- Economical, durable and water resistant
- Combines ultra comfort with high performance and safety
- Withstands even the harshest hospital environments
- Anti-wrinkle material makes it easy to read barcodes and other identifiers
- Durable non-transferable adhesive closure
- High-quality dark & sharp print that doesn't fade
- Easy to print and apply to patients in seconds
- Water-resistant and tamper-evident



WRISTBANDS

INFANT

For babies 0-4 years old provide positive patient identification for newborns using barcode technology to help ensure patient safety. The surface provides best imprint quality, maximum resistance to solvents and seamless printing. Improved direct thermal surface protects print & bar code from water, soap and hand sanitizer. Easy to apply adhesive closure ensures safety and security.

Size of carrier: 29x215 mm, **Size of band:** 19x200 mm, **Information area:** 19x38 mm.

APPLICATIONS

HEALTHCARE

- Patient ID tracking
- Hospital admissions
- Hospital nursery

ENTERTAINMENT & HOSPITALITY

- Amusement parks
- Water parks
- Event admissions
- Night club age verification

OTHER MARKETS

- Crowd control
- Animal ID
- Restricted area access

TECHNICAL INFORMATION

	Infant
Product number	WDI-104-OL
Band width	29 mm
Band length	215 mm
Information area	19mm x 38mm
Printing method	Direct Thermal
Material	Satin
Colour	White
Type of closure	Adhesive
Core size	1,5" / 40 mm
Wound	OUT
Roll O.D.	4,95" / 126mm
Black mark	Yes, at perforation
Bands per roll	200
Packaging	2 rolls/ctn = 400 bands/ctn

PREMIUM WRISTBANDS



Additional features

Phthalate-free, Latex-free, Antimicrobial, Resistant to alcohol, water, soap, hand-sanitizer, and all the most common fluids in the hospital environment.
Certificates: ISO 10993-10 Skin Irritation; Natural Rubber Latex (NRL) free, Bisphenol A (BPA) free, Phthalates (DEHP, DINP, DIDP, DnHP, BBP, DBP, DIBP, DCHP) free.

Wristband safety

Keep wristband clean, dry and don't wear it too tight. When applying a wristband place "one finger under" the band to ensure it is not too tight. When wristbands are worn too tight perspiration, water, soap, lotions etc. caught under the band may cause irritation. If irritation occurs, remove wristband and treat symptoms.

Hand Hygiene: Use proper hand hygiene to reduce the spread of bacteria or other sources of contamination to wristband.

Storage temperature: 22-35°C

Months store: 18 months

Usage: 1-10 days

TECHNICAL SPECIFICATION

